BRAUN

Manual nasal aspirator 1

Helen of Troy

www.BraunHealthcare.com

Certain trademarks used under license from The Procter & Gamble Company or its affiliates.

Kaz Europe Sàrl
Place Chauderon 18
CH-1003 Lausanne
Switzerland

CE

PAP

Made and printed in China

© 2021, All rights reserved.

REF BNA050EU P/N A005537R1 03MAY21

English

Product description

The Manual nasal aspirator is a device that is used to suction mucus from a baby's nasal passages using human suction. The device contains a mouthpiece, tube and mucus receptacle.

Intended use

The Braun Manual nasal aspirator (and associated filters) is intended to be used for providing intermittent suction to remove nasal secretions and mucous in children 0 to 12 years of age, for home use only. The associated filters (in this box and BNF020) can be used for up to 3 months.

Please read all the instructions before using the Braun Manual nasal aspirator, especially this basic safety information.

- The Braun Manual nasal aspirator is designed to clean/extract nasal secretions only. Use only for the intended purpose other than those stated in these instructions for use.
- Not for use of children under 12 years.
- Do not expose the Braun Manual nasal aspirator to direct sunlight.
- This device is designed for domestic use only. Consult a doctor if the child is suffering from a medical infection.
- Do not use other nozzles than those supplied with this device.
- Do not re-use filters.
- Only use Braun Single-use nasal aspirator filters with this device.

CHOKING HAZARD - May contain small detachable parts. Keep out of reach of children under 3 years old.

Cleaning

To clean the Braun Manual nasal aspirator 1, twist off the cup and remove the chamber. Clean it with soap and warm water except the tube. You can also put the detachable part (mouth piece, nose tip and chamber) in the dishwasher except the tube. Alternatively, all parts of the device can be boiled for 15 minutes. This last method is mandatory for the first use.

Replace filter (single use device filter BNF020EU) with every use to maintain hygiene and to prevent cross-contamination.

Guarantee

Please read all instructions before attempting to use this device. Please retain the receipt as proof of date of purchase. The receipt must be presented when making any claim within the relevant guarantee period. Any claim under the guarantee will not be valid without a proof of purchase. Your device is guaranteed for Two Years (2 Years) from date of purchase. This guarantee covers defects in materials or workmanship in this product, or under normal use, defective devices meeting these criteria will be replaced free of charge. The guarantee DOES NOT cover defects or damage resulting from use or failure to follow the user instructions. The guarantee becomes void if the device is opened, tampered with, or used with non-Braun branded parts or accessories, or if used for purposes other than its unauthorized persons. Accessories and consumables are excluded from any warranty. For support requests, please visit www.braunhealthcare.com/uk_en/ or find contact information at the end of this owner's manual.

This Guarantee is applicable to Europe, Russia, Middle East and Africa only.

UK Only: This does not affect your consumer statutory rights.

In Europe, any serious incident that has caused life-threatening injury, surgical intervention) that has occured in relation to this Braun Manual nasal aspirator should be reported to Kaz Europe (see Helpline of Troy Consumer Service Numbers) and to the competent authority of your member state.

REF, LOT and Serial numbers of your product must be provided as they are essential to carry out any follow-up or quality inquiry or claim. The LOT and Serial numbers are located on this sheet.

The production date is given by the LOT number and can be deciphered as explained below:

The first 3 numeric digits in the LOT number represent the day of the year of manufacture.

The next 2 numeric digits represent the last two numbers of the calendar year of manufacture and the letter(s) at the end designate the manufacturer of the product. (E.g.: LOT No.: 1231Zbw this product was made on the day 123, year 2013 at the manufacturer code zbw.)

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Helen of Troy

Creative Department Artwork Specifications

Brand:	Braun
Category:	Nasal Aspirator
Model:	**BNA050EU**
Artwork Part #:	A005537R1
Die Line Part #:	N/A
Agile Desc:	
Subject:	Regulatory Insert Sheet
Region:	EMEA
Flat Size:	MM: **W**750 x **H**320
Folded Size:	MM: **W**75 x **H**106.66
Scale:	1/1
Material:	50# matte text
Page Count:	1
Revision:	1
Date:	**03MAY21**
Release Date:	**05FEB21**
Re-release Date:	**04MAY21**

Colors:

Dielines (Do not print) Spot Colors

Cyan 0%	Magenta 0%	Yellow 0%	Black 100%	PMS XXXX	PMS XXXX

Special Instructions:

Quality Requirement of Artwork and Quality Clarification Process of Artwork Printing: Meet Eng-QS-06&02

Helen of Troy
Creative Services
Marlborough, MA 01752, USA
+1 508 490 7000